

INTELLIGENT  
ULTRASOUND®  
*for smarter scanning*






## HEARTWORKS® | 3D Echo

**Explore complex cardiac anatomy with live 3D volume captures.** Become confident in 3D cardiac image acquisition and manipulation techniques using the world-renowned HeartWorks cardiac simulator - available with all heart models and both TTE and TEE modes.

**3D cardiac image acquisition made simple.** The HeartWorks 3D module guides the student through the steps of data acquisition, offering choices with respect to breath hold and multi-beat acquisition. The simulated datasets allow detailed anatomical inspection of cardiac anatomy and function.

**Learn the principles of 3D dataset manipulation.** 3D datasets can be rotated, zoomed, panned, and cropped with 3 adjustable, colorised orthogonal cropping planes. The simulated datasets allow detailed anatomical inspection of cardiac anatomy and function, as in real life.



**GET IN TOUCH TODAY TO PREBOOK YOUR DEMO**

[hello@intelligentultrasound.com](mailto:hello@intelligentultrasound.com)  
[intelligentultrasound.com](http://intelligentultrasound.com)

